

Ureteral Dilator Sets And Dilators

Application : An access conduit is established and maintained into the bladder through urethra for flexible cystoscopy. 0.038" guide wire is used with this.

Features : Flexible cystoscopy revealed quite non-invasive.


Use : For single use

Size	Product No.	Length
6 FR - 16 FR	NMS-UDS - 6 FR - 16 FR- 70	70 cms.
6 FR - 18 FR	NMS-UDS - 6 FR - 18 FR- 40	40 cms.

Fascial Dilator Sets

Application : It is used for dilation of the ureter , before ureteroscopy or stone manipulation or both. Used with 0.38 " guide wire.

Features : Dialation of ureter with patient safety and comfort.



Ureteral Catheter (Open End)

Application : Used for drainage as well as retrograde pyelogram by open end tip it is placed over a guide wire for negotiation to tortuous ureter.

Features : Smooth passage in tortuous ureter.

Use : For single use

Size	Product No.	Length
3 FR	NMS-UC-3 FR 70 (OE)	70 cms.
4 FR	NMS-UC-4 FR 70(OE)	70 cms.
5 FR	NMS-UC-5 FR 70 (OE)	70 cms.
6 FR	NMS-UC-6 FR 70(OE)	70 cms.

Ureter Catheter (Closed End)

Application : Used for drainage and retrograde pyelogram.


Features : Best in its inherent quality providing drainage and retrograde pyelogram.

Use : For single use

Size	Product No.	Length
3 FR	NMS-UC-3 FR 70 (CE)	70 cms.
4 FR	NMS-UC-4 FR 70 (CE)	70 cms.
5 FR	NMS-UC-5 FR 70 (CE)	70 cms.
6 FR	NMS-UC-6 FR 70 (CE)	70 cms.

NIDHI Meditech Systems

NIDHI Meditech Systems



NIDHI Meditech Systems

(Complete Urological Solutions)

A-207, Tirthraj Complex, Opp. V.S.Hospital,
 B/h. Hasubhai Chambers, Ellisbridge,
 Ahmedabad - 380006. (Gujarat) India.
 Ph. : +91-79-26578322, 26578422, 26553157
 Fax : +91-79-26553157
 Email : sales@nidhimeditech.com
 Internet : www.nidhimeditech.com

NIDHI UROLOGICAL CATHATER



Nidhi Meditech Systems provides genuine endourological products such as ureteral stents. All the stents provided by the company are of the superior quality having the following key features :

1. Durable with superior bio-compatibility
2. Thin wall (provides potency and good drainage.
3. Astraumatic insertion is revealed quite smooth.
4. Minimized tissue irritation with easy but rapid astraumatic insertion.
5. Cost effective.



Double Pigtail Ureteral Stents

(Both end Open as well as One end Close)

Application : Applied for internal drainage from the ureteropelvic junction to the bladder, as a provisional measure.

Features : Made of soft telecoflex polyurethane for smooth passage radio-opaque implant.

Use : For single use

CE = One End Close, BEO = Both End Open

Size	Product No.	Length	Guide Wire..
4.5 FR	NMS-DJ-4.5 FR 24 (CE)	24 cms.	0.025"
4.5 FR	NMS-DJ-4.5 FR 24 (BEO)	24 cms.	0.025"
4.5 FR	NMS-DJ-4.5 FR 26 (CE)	26 cms.	0.025"
4.5 FR	NMS-DJ-4.5 FR 26 (BEO)	26 cms.	0.025"
4.8 FR	NMS-DJ-4.8 FR 24 (CE)	24 cms.	0.032"
4.8 FR	NMS-DJ-4.8 FR 24 (BEO)	24 cms.	0.032"
4.8 FR	NMS-DJ-4.8 FR 26 (CE)	26 cms.	0.032"
4.8 FR	NMS-DJ-4.8 FR 26 (BEO)	26 cms.	0.032"
5 FR	NMS-DJ- 5 FR 24 (CE)	24 cms.	0.035"
5 FR	NMS-DJ-5 FR 24 (BEO)	24 cms.	0.035"
5 FR	NMS-DJ- 5 FR 26 (CE)	26 cms.	0.035"
5 FR	NMS-DJ-5 FR 26 (BEO)	26 cms.	0.035"
5 FR	NMS-DJ- 5 FR 28 (CE)	28 cms.	0.035"
5 FR	NMS-DJ-5 FR 28 (BEO)	28 cms.	0.035"
5.5 FR	NMS-DJ- 5.5 FR 26 (CE)	26 cms.	0.035"
5.5 FR	NMS-DJ-5.5 FR 26 (BEO)	26 cms.	0.035"
5.5 FR	NMS-DJ- 5.5 FR 28 (CE)	28 cms.	0.035"
5.5 FR	NMS-DJ-5.5 FR 28 (BEO)	28 cms.	0.035"
6 FR	NMS-DJ- 6 FR 26 (CE)	26 cms.	0.038"
6 FR	NMS-DJ- 6 FR 26 (BEO)	26 cms.	0.038"
6 FR	NMS-DJ- 6FR - 28 (CE)	28 cms.	0.038"
6 FR	NMS-DJ- 6FR - 28 (BEO)	28 cms.	0.038"


Stone Basket - Helical (Dormia Basket)

Application : To remove multiple stone fragments from ureter.

Features : Maximum patient safety and compatibility. Easy to operate.

Use : For single use

Size (In FR)	Product No.	Length
3 / 4 WIRE	NMS-DB-H-3FR-4FR-70	70 cms.
4 / 4 WIRE	NMS-DB-H-4FR-4FR-70	70 cms.
4 / 5 WIRE	NMS-DB-H-4FR-5FR-70	70 cms.
5 / 4 WIRE	NMS-DB-H-5FR-4FR-70	70 cms.
5 / 5 WIRE	NMS-DB-H-5FR-5FR-70	70 cms.



Flatwire Stone Retrieval Basket

Application : It is used for removing the large stones for Ureter.

Features : Retrieving the large stones, smoothly.

Use : For single use

Percutaneous Pigtail Nephrostomy Catheter With Needle

(PCN Catheter With Needle)

Application : By percutaneous placement of pigtail catheter, bladder drainage is provided.

Features : Easy outlet for better drainage.

Use : For single use



Size	Product No.	Length
6 FR	NMS-PCNC-N-6 FR- 22	22 cms.
6 FR	NMS-PCNC-N-6 FR 30	30 cms.
8 FR	NMS-PCNC-N-8 FR 22	22 cms.
8 FR	NMS-PCNC-N-8 FR - 30	30 cms.
9 FR	NMS-PCNC-N-9 FR 22	22 cms.
9 FR	NMS-PCNC-N-9 FR 30	30 cms.
10 FR	NMS-PCNC-N-10FR 22	22 cms.
10 FR	NMS-PCNC-N-10FR 30	30 cms.
12 FR	NMS-PCNC-N-12FR 22	22 cms.
12 FR	NMS-PCNC-N-12FR 30	30 cms.
14 FR	NMS-PCNC-N-14FR 22	22 cms.
14 FR	NMS-PCNC-N-14FR 30	30 cms.

Malecot Nephrostomy Catheter (W/O Needle)

Application : It is used to provide nephrostomy drainage by percutaneous placement of malecot catheter in the renal pelvis.

Features : Sterile in peel open package along with Urine bag connector.

Use : For single use



Size	Product No.	Length
8 FR	NMS-MLC-N-8 FR-30	30 cms.
8 FR	NMS-MLC-8 FR-30 (Without Needle)	30 cms.
10 FR	NMS-MLC-N-10FR 30	30 cms.
10 FR	NMS-MLC-10FR 30 (Without Needle)	30 cms.
12 FR	NMS-MLC-N-12 FR 30	30 cms.
12 FR	NMS-MLC-12 FR 30 (Without Needle)	30 cms.
14 FR	NMS-MLC-N-14 FR - 30	30 cms.
14 FR	NMS-MLC-14 FR - 30 (Without Needle)	30 cms.
16 FR	NMS-MLC-N-16 FR 30	30 cms.
16 FR	NMS-MLC-16 FR 30 (Without Needle)	30 cms.
18 FR	NMS-MLC-N-18 FR 30	30 cms.
18 FR	NMS-MLC-18 FR 30 (Without Needle)	30 cms.

Percutaneous Suprapubic Catheter Sets (Malecot Catheter With Needle)

Application : It is used to provide bladder drainage by percutaneous placement of malecot catheter.

Features : Bladder drainage with full patient comfort.

Use : For single use

Pigtail Nephrostomy Catheter (PCN Catheter) Without Needle

Application : By direct puncture it provides temporary (even permanent) drainage of urine, from the kidney. This is done in presence of distal obstruction.

Features : Provides drainage from kidney, smoothly.

Use : For single use

Size	Product No.	Length
6 FR	NMS-PCNC-N-6 FR- 22	22 cms.
6 FR	NMS-PCNC-N-6 FR 30	30 cms.
8 FR	NMS-PCNC-N-8 FR 22	22 cms.
8 FR	NMS-PCNC-N-8 FR - 30	30 cms.
9 FR	NMS-PCNC-N-9 FR 22	22 cms.
9 FR	NMS-PCNC-N-9 FR 30	30 cms.
10 FR	NMS-PCNC-N-10FR 22	22 cms.
10 FR	NMS-PCNC-N-10FR 30	30 cms.
12 FR	NMS-PCNC-N-12FR 22	22 cms.
12 FR	NMS-PCNC-N-12FR 30	30 cms.
14 FR	NMS-PCNC-N-14FR 22	22 cms.
14 FR	NMS-PCNC-N-14FR 30	30 cms.